

Bierley

MonoMouse-RM™



The MonoMouse-RM provides three modes for reading, namely regular mode, reverse mode (pure white text on black background), and enhanced mode (pure black text on white background). With quick installation by plugging the MonoMouse-RM into any regular television and pressing the large blue button, the magnifier is ready to be used and your television becomes a powerful low vision aid. Operation couldn't be easier because you only need to switch to the video channel using your remote control and when you have finished reading, switch back to your favourite channel.

You can leave your MonoMouse-RM connected to the television so that it is always ready to use.

**Non-Glare
Image**

The ColorMouse-RM offers all the features of our TV range in one product! Full colour for magazine, greyscale for books and Reverse Mode for Macular Degeneration.

**White on
Black**

**Black on
White**

Full Colour

**1
Year
Worldwide
Warranty**

www.bierley.com

- Compatible with any TV using video input
- High contrast image, less glare
- Ergonomic and easy to use
- High contrast negative image (pure white text on black background)
- High contrast positive image (pure black text on white background)
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- RCA connection
- 6 feet of easy store cable
- Durable, CE approved, built to last

**CALL NOW
03 9877 7956**