

# Bierley



## MonoMouse-RM™



The MonoMouse-RM provides three modes for reading, namely regular mode, reverse mode (pure white text on black background), and enhanced mode (pure black text on white background). With quick installation by plugging the MonoMouse-RM into any regular television and pressing the large blue button, the magnifier is ready to be used and your television becomes a powerful low vision aid. Operation couldn't be easier because you only need to switch to the video channel using your remote control and when you have finished reading, switch back to your favorite channel.

You can leave your MonoMouse-RM connected to the television so that it is always ready to use.

**Non-Glare  
Image**

**Reverse Mode makes the MonoMouse-RM the choice for anyone with Macular Degeneration. 100% contrast provides the highest level of clarity possible.**

**White on  
Black**

**Black on  
White**

- Compatible with any TV using video input
- High contrast image, less glare
- Ergonomic and easy to use
- High contrast negative image (pure white text on black background)
- High contrast positive image (pure black text on white background)
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- RCA connection
- 6 feet of easy store cable
- Durable, CE approved, built to last

**1  
Year  
Worldwide  
Warranty**

[www.bierley.com](http://www.bierley.com)

**CALL NOW**

**1 411 6486-2500**