

# Bierley

## MonoMouse-RM™



**Non-Glare  
Image**

The MonoMouse-RM provides three modes for reading, namely regular mode, reverse mode (pure white text on black background), and enhanced mode (pure black text on white background). With quick installation by plugging the MonoMouse-RM into any regular television and pressing the large blue button, the magnifier is ready to be used and your television becomes a powerful low vision aid. Operation couldn't be easier because your television will automatically switch to the correct channel for using your MonoMouse-RM and will automatically switch back to the channel that you were watching once you switch the MonoMouse-RM off. This feature is unique to our European SCART products.

You can leave your MonoMouse-RM connected to the television so that it is always ready to use.

**Reverse Mode makes the MonoMouse-RM the choice for anyone with Macular Degeneration. 100% contrast provides the highest level of clarity possible.**

**White on  
Black**

**Black on  
White**

- Compatible with any TV using video input
- High contrast image, less glare
- Ergonomic and easy to use
- High contrast negative image (pure white text on black background)
- High contrast positive image (pure black text on white background)
- 24x Magnification on a 20" TV :(High)
- 13x Magnification on a 20" TV :(Low)
- SCART connection with Auto-Detect
- 9 feet of easy store cable
- Durable, CE approved, built to last

**1  
Year  
Worldwide  
Warranty**

[www.bierley.com](http://www.bierley.com)

**CALL NOW  
0800 0430282**